

SPLINKERETTE Protocol

Steps 1-3 are currently conducted in 96-well microtiter plates whereas steps 4-6 are carried out in 384-well format.

1. Genomic DNA (gDNA) preparation from frozen ES cell clone pellets

Lysis buffer: 10mM Tris/HCl pH 7.5
 10mM EDTA
 10mM NaCl
 0.5% SDS
 1mg/ml Proteinase K (to be added freshly)

- 1.1. starting material: A 96-well microtiter plate containing PBS-washed frozen ES cell clone pellets
- 1.2. add 40 µl lysis buffer to each well
- 1.3. seal 96-well plate tightly (e.g. using Microseal A Film [Biozym, Art.Nr. 621715])
- 1.4. incubate at 60°C for 2h
- 1.5. cool down 96-well plate to room temperature, spin-down condensed liquid (20°C, 2min, 4000g)
- 1.6. add 80 µl icecold 75mM (w/v) NaCl/EtOH suspension
- 1.7. incubate 1 h at room temperature for gDNA precipitation
- 1.8. slowly invert plate over a sink to discard the ethanolic solution. The precipitated nucleic acid should remain attached to the base of the wells and may be visible against bright light.
- 1.9. allow the remaining ethanol to drain from the plate on a bed of paper towels
- 1.10. wash gDNA pellet 3x with 150 µl 70% EtOH each, being careful not to dislodge the precipitate
- 1.11. dry gDNA pellet in speed-vac for 5-10min until the last traces of ethanol have evaporated
- 1.12. elute genomic DNA with 40 µl 5mM Tris/HCl pH 8.0
- 1.13. dissolve gDNA over night (12-16 h) at room temperature

2. Restriction digest and ligation

2.1. prepare restriction master mix (final reaction volume: 20 μ l)

2.1.1. for FlipROSAbetageo vectors (V08, V09, V10, V11, V12, V13, V14):

	for 1 sample
10x NEB buffer 2	2.0 μ l
<i>Bst</i> YI (NEB, 10U/ μ l)	0.5 μ l
Bidest-water	8.5 μ l

2.1.2. for ROSAbetageo vectors (V04, V16, V17):

	for 1 sample
10x NEB buffer 3	2.0 μ l
100x BSA (NEB)	0.2 μ l
<i>Apo</i> I (NEB, 10U/ μ l)	0.5 μ l
Bidest-water	8.3 μ l

2.2. pipet 11.0 μ l restriction mix to each sample

2.3. transfer 9.0 μ l of gDNA to each sample

2.4. incubate:

2.5 h at 60°C

20 min at 80°C

cool down to 4°C

2.5. prepare Splinkerette: for 1 sample

Primer SpAa, 10 μ M: 0.15 μ l

Primer SpBb, 10 μ M: 0.15 μ l

SuRE buffer M (Roche, Basel): 0.05 μ l

Bidest water: 0.65 μ l

2.6. incubate Splinkerette mixture for 4 min at 97°C

slowly let cool down on benchtop for annealing

2.7. prepare ligation master mix as follows (final reaction volume: 30 µl)

	for 1 sample
Splinkerette:	1.0 µl
10x DNA Ligation buffer (NEB)	3.0 µl
T4-DNA Ligase (NEB, 400U/µl)	1.0 µl
Bidest-water	5.0 µl

2.8. pipette 10.0 µl ligation mix to each sample

2.9. incubate over night (12-16h) at 16°C

3. Ligated genomic DNA purification using Millipore-Multiscreen-PCR-Filter-plates

- 3.1. add 25 µl Bidest-water to Multiscreen-PCR-plate to wet the membrane
- 3.2. dilute the 30 µl gDNA with 150 µl Bidest-water and transfer it into the Multiscreen-PCR-plate
- 3.3. place the Multiscreen PCR-plate on top of the Multiscreen-manifold, place a plate for waste
- 3.4. apply vakuüm at 15 inches for 10 min
- 3.5. add 40 µl 5mM Tris pH 8.0
- 3.6. retrieve clean ligated DNA by pipetting up and down several times
- 3.7. transfer cleaned DNA into a new 96-well plate
- 3.8. merge cleaned DNA of 4 x 96-well plates to 1 x 384-well plate
by transferring 5 µl of clean DNA according to the following pipetting scheme:

- 384 well A01 corresponds to 1st 96 well A01
- 384 well A02 corresponds to 2nd 96 well A01
- 384 well B01 corresponds to 3rd 96 well A01
- 384 well B02 corresponds to 4th 96 well A01
- ...

3.9. continue to first round PCR

	P	O	N	M	L	K	J	I	H	G	F	E	D	C	B	A
1																
2																
3																
4																
5																
6																
7																
8																
9																
10																
11																
12																
13																
14																
15																
16																
17																
18																
19																
20																
21																
22																
23																
24																

4. First round PCR

primer pair for 5' side: Sp0F/5'Sp0R

primer pair for 3' side: Sp0F/3'Sp0R

<u>component</u>	<u>amount (µl)</u>
gDNA	5.00
sterilized, distilled water	6.17
10x reaction buffer (incl. 1.5mM Mg ²⁺)	1.50
25mM MgCl ₂	0.15
2mM dNTP mix	1.50
primer Sp0F (10pmol/µl)	0.30
primer Sp0R (10pmol/µl)	0.30
<u>Eppendorf Taq Polymerase (0.04 U/µl)</u>	<u>0.08</u>
final volume	15.00

program: 94°C/75s + 2x(94°C/20s; 64°C/15s) + 30x(94°C/20s; 58-64°C/15s; 72°C/120s) + 72°C/7mins

5. Second round PCR

primer pair for 5' side: Sp1F/5'Sp1R

primer pair for 3' side: Sp1F/3'Sp1R

<u>component</u>	<u>amount (µl)</u>
first round PCR product	1.00 (or minute amount by stamp)
sterilized, distilled water	6.42
10x reaction buffer	1.00
25mM MgCl ₂	0.10
2 mM dNTP mix	1.00
primer Sp1F (10pmol/µl)	0.20
primer Sp1R (10pmol/µl)	0.20
<u>Taq Polymerase (Eppendorf)</u>	<u>0.08</u>
final volume	10.00

program: 94°C/60s + 30x(94°C/20s; 60°C/15s; 72°C/60s) + 72°C/7mins

analyze 2-3µl of secondary amplification product on 1.4% agarose gel or proceed directly to sequencing

6. BigDye Terminator Sequencing PCR

<u>component</u>	<u>amount (µl)</u>
second round PCR product	1.00
sterilized, distilled water	7.00
seq. primer (10pmol/µl)	1.00
<u>BigDye Terminator Mix</u>	<u>1.00</u>
final volume	10.00

program: 96°C/60s + 25x(96°C/10s; 55°C/5s; 60°C/240s)

7. Primers

Number	Sequence (5' -> 3')	Primer type	Tm value
7.1. ROSAbetageo vectors (V04, V16, V17)			
40	CGAAGAGTAACCGTTGCTAGGAGAGACCGTGGCTG AATGAGACTGGTGTGCGACACTAGTGG	Splinkerette SpAa	
62	AATCCACTAGTGTGCGACACCAGTCTCTAATTTTTTTT TTTTCAAAAAA	Splinkerette SpBb	
38	CGAAGAGTAACCGTTGCTAGGAGAGACC	Sp0F	62°C
56	AGGAAACCCTGGACTACTGC	5' Sp0R	62°C
58	GACACAGATAAGTTGCTGGC	3' Sp0R	60°C
39	GTGGCTGAATGAGACTGGTGTGCGAC	Sp1F	61°C
47	GATGGAACAGCTAGAGAACC	5' Sp1R	60°C
43	CTCAGTTATGTATTTTTCCATGCCTTGC	3' Sp1R	59°C
21	TTGTGGTCTCGCTGTTCCCTTGGG	5' Seq	62°C
37	GCTAGCTTGCCAAACCTACAGGTGG	3' Seq	61°C
7.2. FlipROSAbetageo vectors (V08, V09, V10, V11, V12, V13, V14)			
40	CGAAGAGTAACCGTTGCTAGGAGAGACCGTGGCTG AATGAGACTGGTGTGCGACACTAGTGG	Splinkerette SpAa	
41	GATCCCACTAGTGTGCGACACCAGTCTCTAATTTTTTTT TTTTCAAAAAA	Splinkerette SpBb	
38	CGAAGAGTAACCGTTGCTAGGAGAGACC	Sp0F	62°C
33	CGACCAGCTGTGCGCATAGTG	5' Sp0R	59°C
35	AGTCATAGACACTAGACAATCGG	3' Sp0R	60°C
39	GTGGCTGAATGAGACTGGTGTGCGAC	Sp1F	61°C
34	TTTGGCAAGCTAGCACAACC	5' Sp1R	60°C
36	CAGTCAATCGGAGGACTGGCG	3' Sp1R	59°C
21	TTGTGGTCTCGCTGTTCCCTTGGG	5' Seq	62°C
37	GCTAGCTTGCCAAACCTACAGGTGG	3' Seq	61°C