

## FRT-PCR Protocol

1. Starting material: A 96-well microtiter plate containing ~ 30  $\mu$ l genomic ES cell DNA.
2. Add 40  $\mu$ l 5 mM Tris/HCl pH 8.0 to each well.
3. Merge 1:2 diluted DNA of 4 x 96-well plates to 1 x 384-well plate by transferring 5  $\mu$ l of diluted DNA according to the following pipetting scheme:

384 well A01 corresponds to 1<sup>st</sup> 96 well A01

384 well A02 corresponds to 2<sup>nd</sup> 96 well A01

384 well B01 corresponds to 3<sup>rd</sup> 96 well A01

384 well B02 corresponds to 4<sup>th</sup> 96 well A01

...

4. Continue to FRT-PCR:

Component	Amount ( $\mu$ l)
Genomic DNA	5,00
Sterilised, distilled water	2,12
10x reaction buffer (incl. 1,5 mM Mg <sup>2+</sup> )	1,00
2 mM dNTP mix	1,00
Sense primer 2 (10 pmol/ $\mu$ l)	0,40
Antisense primer 4 (10 pmol/ $\mu$ l)	0,40
Eppendorf Taq Polymerase (5U/ $\mu$ l)	0,08
Final volume	10,00

Add 5  $\mu$ l of PCR mix to each 384 well.

Program: 94°C/120sec + 34x(96°C/15sec; 55°C/15sec; 72°C/60sec) + 72°C/7min

5. Analyse 2  $\mu$ l of amplification product on 1,4% agarose gel or proceed directly to sequencing.

## 6. BigDye (v3) Terminator Sequencing PCR:

Component	Amount ( $\mu$ l)
FRT-PCR product	1,00
Sterilised, distilled water	7,00
Seq. primer 1 (10 pmol/ $\mu$ l)	1,00
Big Dye Terminator Mix	1,00
Final volume	10,00

Program: 96°C/60sec + 25x(96°C/10sec; 56°C/5sec; 70°C/240sec)

## 7. Primers

Number	Sequence (5'→ 3')	Primer type	Tm value
2	CGAGTGAGGGGTTGTGCTAG	FRT Fwd	61°C
4	CCTATAGTGAGTCGTATTCTCCC G	FRT Rev	63°C
1	CACTATGCGCACAGCTGGTCCG	Seq. Fwd	62°C